

---

# Anamorphic prism pair (White 440-640 nm)

Optical parts and Scanner sets



---

Anamorphic prism pair is used to expand or reduce the beam size in one direction only. This pair is optimized for all visible light from 440 nm to 640 nm. The magnification can be well tuned continuously by rotating the prisms, although the coating is optimised for magnification 1:2. There are 2 dimensions available: height 16 mm, width 12 mm, thickness near 8 mm or height 16mm, width 10mm, thickness near 7mm. Please specify the preferred size in you order.